


Xtrac EC

Aspiration Catheter

SIS  MEDICAL

Swiss | Interventional | Systems



- a) Marker band 3 mm from distal tip
- b) Aspiration lumen (elliptical shape)
- c) 0.014" guidewire
- d) Rapid exchange part, 25 cm, highly flexible and hydrophilic coated
- e) Stiffening stylet
- f) Rapid exchange wire exit

SIS MEDICAL

Swiss | Interventional | Systems

Hungerbühlstrasse 12a | CH-8500 Frauenfeld
 Tel. +41 (0)52 245 09 90 | Fax +41 (0)52 245 09 91
 info@sis-medical.com | www.sis-medical.com



Key Features:

- Easy and fast placement provided by a long, hydrophilic coated rapid exchange part (25 cm) and a supporting stylet which is significantly improving the pushability and kink resistance of the catheter.
- Precise positioning supported by very good visibility of tip marker 3 mm from the distal, atraumatic catheter end.
- Large aspiration lumen providing a high aspiration capacity.
- Large and elliptical aspiration hole shape facilitating the entry of large pieces of thrombus.
- 60 ml syringe saving time and increasing aspiration effectiveness.

Technical Specifications:

6F version compatible with Guding Catheter 6F with an inner lumen of min 0.071"

7F version compatible with Guding Catheter 7F with an inner lumen of min 0.081"

• Method	distal tip aspiration
• Total length	145.00 cm
• Rx length	25.00 cm
• Hydrophilic coating on the distal 25 cm of the catheter	
• Stiffening stylet with rounded distal edge	
• Length of the stiffening stylet	126.50 cm
• Tip length distal to radiopaque marker	3.00 mm
• Distal shaft diameter	6F = 1.70 mm / 7F = 1.98 mm
• Proximal shaft diameter	6F = 1.30 mm / 7F = 1.60 mm
• Aspiration capacity (water)	6F = 60 ml/45 sec or 1.33 ml/sec 7F = 60 ml/22 sec or 2.77 ml/sec
• Min. inner diameter	6F = 0.041" / 7F = 0.053"

The Aspiration Catheter Set:

- 1x Aspiration Catheter 6F or 7F
- 1x Stiffening Stylet
- 1x Aspiration Syringe 60 ml
- 2x Filter Baskets
- 1x Line Extension with one way stopcock
- 1x Instruction for Use

Ordering Information:

1. Aspiration Catheter 6F ref nr.: 04030271
2. Aspiration Catheter 7F ref nr.: 04030272



SIS MEDICAL

Swiss | Interventional | Systems

Hungerbühlstrasse 12a | CH-8500 Frauenfeld
Tel. +41 (0)52 245 09 90 | Fax +41 (0)52 245 09 91
info@sis-medical.com | www.sis-medical.com

Copyright © by SIS Medical AG.
Currently not available for sale in USA.
Subject to alterations due to technical advances.
All rights reserved.