

NIC 1.1 hydro



1.1 MM LOW-PROFILE
NON-COMPLIANT PTCA
DILATATION CATHETER

21 atm

High RBP

Tip

Unique tip design to facilitate navigation

Lubricious hydrophilic coating

Marker

One Pt/Ir central marker

LOW PROFILE HYDROPHILIC COATING



SWISS MADE


www.sis-medical.com

Manufactured by SIS Medical AG.
All rights reserved.

NIC 1.1 hydro

LOW PROFILE NON-COMPLIANT PTCA DILATATION CATHETER

Balloon Diameter, mm	Balloon Length, mm	Art. Number	UDI-DI-Number (GTIN)
1.1	6	110-006-134	7640142573098
1.1	10	110-010-134	7640142573104
1.1	15	110-015-134	7640142573128
1.1	20	110-020-134	7640142573142



TECHNICAL SPECIFICATIONS

Nominal Pressure	10 atm
Rated Burst Pressure (RBP)*	21 atm
Guiding Catheter Compatibility	5F
Guide Wire Compatibility	0.014" (0.36 mm)
Marker position	Single Pt/Ir central marker
Catheter Usable Length	140 cm
Catheter Coating	Hydrophilic

*Do not exceed RBP

NIC 1.1 hydro is not approved for sale in the USA.

SIS MEDICAL

Hungerbühlstrasse 12a
CH-8500 Frauenfeld
Switzerland
+41 52 245 09 90
sales@sis-medical.com
www.sis-medical.com

GERMANY
Distributed by: SIS Medical GmbH
Kirchstrasse 36, 72127, Kusterdingen, Germany
+49 7127-799 8830
orders-de@sis-medical.com

Illustrations are not to be recognized as engineering drawings.

CE marked product.

LA002639-001 03/2025 NIC 1.1 hydro brochure