


Xtrac

**ASPIRATION PERFORMANCE SPECIFICALLY TAILORED FOR
THROMBUS REMOVAL**



ASPIRATION CATHETER

Xtrac

INTENDED USE

Indicated for use in the central and peripheral circulatory system, including saphenous vein grafts to contain and aspirate embolic material (thrombus/debris) while performing percutaneous transluminal coronary angioplasty, percutaneous transluminal angioplasty, and/or stenting procedures.

Reference Description Hydrophilic coating.

LA 002651-001 Brochure Xtrac

Art. Numbers	Description	Hydrophilic coating	Stylet
04030271	Xtrac EC 6F Stiff Sterile	yes	yes
04030272	Xtrac EC 7F Stiff Sterile	yes	yes

Important information: Prior to use, refer to the “instructions for use” supplied with these devices for indications, contraindications, side effects, suggested procedures, warnings and precautions.

For healthcare professionals only. Class III medical device.



Package Content:



- > Thermoformed support
- > Xtrac catheter and stylet
- > Syringe Vaclock
- > Cell stainer
- > Extension line

SIS MEDICAL



Hungerbühlstrasse 12a
CH-8500 Frauenfeld
Switzerland
+41 52 245 09 90
sales@sis-medical.com
www.sis-medical.com

Manufactured by Arthesys
4 rue Rene Razel
91400 Saclay
France



All rights reserved.

Illustrations are not to be recognized as engineering drawings.