

# INFLATION DEVICE



**SUPPORTS ULTRA-HIGH PRESSURE  
INFLATION UP TO 55 atm**

**QUICK LATCH DESIGN**

**40 atm**

**55 atm**

Designed for balloon catheter dilation, including vertebral cavity creation and coronary interventions. It features a stable, ergonomic handle with a quick latch mechanism allowing rapid deflation even at maximum pressure.

Art. Numbers	Description	atm	Type
96463	Inflation Device	55	Quick Latch
96346	Inflation Device	40	Quick Latch



**SIS MEDICAL**  
 Hungerbühlstrasse 12a  
 CH-8500 Frauenfeld  
 Switzerland  
 +41 52 245 09 90  
 sales@sis-medical.com  
 www.sis-medical.com

Manufactured by  
 Atrion Medical Products, Inc.  
 1426 Curt Francis Road, Arab, AL 35016  
 www.atrionmedical.com  
 All rights reserved.



Illustrations are not to be recognized as engineering drawings.